

Rat Ccdc82/Coiled-coil domain-containing protein 82 ELISA Kit

Catalog No: #EK20055

Package Size: #EK20055 1x96T

Description

Product Name	Rat Ccdc82/Coiled-coil domain-containing protein 82 ELISA Kit
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Coiled-coil domain-containing protein 82
Target Name	Coiled-coil domain-containing protein 82
Accession No.	Q66H73

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.182ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Ccdc82

Note: This product is for in vitro research use only and is not intended for use in humans or animals.