

Rat Crry/Complement regulatory protein Crry ELISA Kit

Catalog No: #EK20056

Package Size: #EK20056 1x96T

Description

Product Name	Rat Crry/Complement regulatory protein Crry ELISA Kit
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Complement regulatory protein Crry
Target Name	Complement regulatory protein Crry
Other Names	Antigen 5I2
Accession No.	Q63135

Application Details

Assay Type: Sandwich Detect Range: 0.312-20 ng/mL Sensitivity: 0.183 ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100 µl Assay Time: 3-5 h Detection wavelength: 450 nm gene_full_name: Crry

Note: This product is for in vitro research use only and is not intended for use in humans or animals.