

Human ASGR2/Asialoglycoprotein receptor 2 ELISA Kit

Catalog No: #EK20057

Package Size: #EK20057 1x96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ASGR2/Asialoglycoprotein receptor 2 ELISA Kit
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Asialoglycoprotein receptor 2
Target Name	Asialoglycoprotein receptor 2
Other Names	ASGP-R 2,C-type lectin domain family 4 member H2,Hepatic lectin H2,HL-2,CLEC4H2
Accession No.	P07307

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:ASGR2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.