

Human BMP5/Bone morphogenetic protein 5 ELISA Kit

Catalog No: #EK20058

Package Size: #EK20058 1x96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human BMP5/Bone morphogenetic protein 5 ELISA Kit
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Bone morphogenetic protein 5
Target Name	Bone morphogenetic protein 5
Other Names	BMP-5
Accession No.	P22003

Application Details

Assay Type: Sandwich Detect Range: 0.312-20 ng/mL Sensitivity: 0.156 ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100 µL Assay Time: 3-5 h Detection wavelength: 450 nm gene_full_name: BMP5

Note: This product is for in vitro research use only and is not intended for use in humans or animals.