

Human COL4A4/Collagen alpha-4(IV) chain ELISA Kit

Catalog No: #EK20062

Package Size: #EK20062 1x96T

Description

Product Name	Human COL4A4/Collagen alpha-4(IV) chain ELISA Kit
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Collagen alpha-4(IV) chain
Target Name	Collagen alpha-4(IV) chain
Accession No.	P53420

Application Details

Assay Type: Sandwich Detect Range: 0.312-20 ng/mL Sensitivity: 0.156 ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100 µl Assay Time: 3-5 h Detection wavelength: 450 nm gene_full_name: COL4A4

Note: This product is for in vitro research use only and is not intended for use in humans or animals.