

Mouse Adc/Antizyme inhibitor 2 ELISA Kit

Catalog No: #EK20063

Package Size: #EK20063 1x96T

Description

Product Name	Mouse Adc/Antizyme inhibitor 2 ELISA Kit
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Antizyme inhibitor 2
Target Name	Antizyme inhibitor 2
Other Names	AZI2, Arginine decarboxylase-like protein, ADC, ARGDC
Accession No.	Q8BVM4

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.068ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Adc

Note: This product is for in vitro research use only and is not intended for use in humans or animals.