

Rat Axin2/Axin-2 ELISA Kit

Catalog No: #EK20065

Package Size: #EK20065 1x96T



Description

Product Name	Rat Axin2/Axin-2 ELISA Kit
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Axin-2
Target Name	Axin-2
Other Names	Axin-like protein,Axil,Axis inhibition protein 2,Conductin
Accession No.	O70240

Application Details

Assay Type: SandwichDetect Range:0.625-40ng/mL Sensitivity:0.312ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Axin2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.