

Rat Afm/Afamin ELISA Kit

Catalog No: #EK20069

Package Size: #EK20069 1x96T

Description

Product Name	Rat Afm/Afamin ELISA Kit
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Afamin
Target Name	Afamin
Other Names	Alpha-albumin, Alpha-Alb
Accession No.	P36953

Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.38ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Afm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.