

Human ANTXR1/Anthrax toxin receptor 1 ELISA Kit

Catalog No: #EK20072

Package Size: #EK20072 1x96T



Description

Product Name	Human ANTXR1/Anthrax toxin receptor 1 ELISA Kit
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Anthrax toxin receptor 1
Target Name	Anthrax toxin receptor 1
Other Names	Tumor endothelial marker 8,ATR,TEM8
Accession No.	Q9H6X2

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.082ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:ANTXR1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.