

Mouse Ass1/Argininosuccinate synthase ELISA Kit

Catalog No: #EK20075

Package Size: #EK20075 1x96T



Description

Product Name	Mouse Ass1/Argininosuccinate synthase ELISA Kit
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Argininosuccinate synthase
Target Name	Argininosuccinate synthase
Other Names	Citrulline--aspartate ligase,Ass
Accession No.	P16460

Application Details

Assay Type: SandwichDetect Range:78-5000pg/mL Sensitivity:39.4pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Ass1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.