

Human ABCC11/ATP-binding cassette sub-family C member 11 ELISA Kit

Catalog No: #EK20077

Package Size: #EK20077 1x96T

Description

Product Name	Human ABCC11/ATP-binding cassette sub-family C member 11 ELISA Kit
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human ATP-binding cassette sub-family C member 11
Target Name	ATP-binding cassette sub-family C member 11
Other Names	Multidrug resistance-associated protein 8,MRP8
Accession No.	Q96J66

Application Details

Assay Type: SandwichDetect Range:78-5000pg/mL Sensitivity:39.6pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:ABCC11

Note: This product is for in vitro research use only and is not intended for use in humans or animals.