

# Human BCL10/B-cell lymphoma/leukemia 10 ELISA Kit

Catalog No: #EK20078

Package Size: #EK20078 1x96T

## Description

|                    |                                                                                                                                                                                                                                                                                               |
|--------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name       | Human BCL10/B-cell lymphoma/leukemia 10 ELISA Kit                                                                                                                                                                                                                                             |
| Applications       | ELISA                                                                                                                                                                                                                                                                                         |
| Species Reactivity | Hu                                                                                                                                                                                                                                                                                            |
| Specificity        | Natural and recombinant Human B-cell lymphoma/leukemia 10                                                                                                                                                                                                                                     |
| Target Name        | B-cell lymphoma/leukemia 10                                                                                                                                                                                                                                                                   |
| Other Names        | B-cell CLL/lymphoma 10,Bcl-10,CARD-containing molecule enhancing NF-kappa-B,CARD-like apoptotic protein,hCLAP,CED-3/ICH-1 prodomain homologous E10-like regulator,CIPER,Cellular homolog of vCARMEN,cCARMEN,Cellular-E10,c-E10,Mammalian CARD-containing adapter molecule E10,mE10,CIPER,CLAP |
| Accession No.      | O95999                                                                                                                                                                                                                                                                                        |

## Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.086ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene\_full\_name:BCL10

Note: This product is for in vitro research use only and is not intended for use in humans or animals.