

Mouse Bmf/Bcl-2-modifying factor ELISA Kit

Catalog No: #EK20080

Package Size: #EK20080 1x96T

Description

Product Name	Mouse Bmf/Bcl-2-modifying factor ELISA Kit
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Bcl-2-modifying factor
Target Name	Bcl-2-modifying factor
Accession No.	Q91ZE9

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.087ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Bmf

Note: This product is for in vitro research use only and is not intended for use in humans or animals.