

Rat Bmf/Bcl-2-modifying factor ELISA Kit

Catalog No: #EK20081

Package Size: #EK20081 1x96T

Description

Product Name	Rat Bmf/Bcl-2-modifying factor ELISA Kit
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Bcl-2-modifying factor
Target Name	Bcl-2-modifying factor
Accession No.	Q8K589

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.175ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Bmf

Note: This product is for in vitro research use only and is not intended for use in humans or animals.