

Human BRMS1/Breast cancer metastasis-suppressor 1 ELISA Kit

Catalog No: #EK20082

Package Size: #EK20082 1x96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human BRMS1/Breast cancer metastasis-suppressor 1 ELISA Kit
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Breast cancer metastasis-suppressor 1
Target Name	Breast cancer metastasis-suppressor 1
Accession No.	Q9HCU9

Application Details

Assay Type: Sandwich Detect Range: 78-5000 pg/mL Sensitivity: 39.7 pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: BRMS1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.