

Human CASP6/Caspase-6 ELISA Kit

Catalog No: #EK20085

Package Size: #EK20085 1x96T



Description

Product Name	Human CASP6/Caspase-6 ELISA Kit
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Caspase-6
Target Name	Caspase-6
Other Names	CASP-6,Apoptotic protease Mch-2,MCH2
Accession No.	P55212

Application Details

Assay Type: SandwichDetect Range:31.2-2000pg/mL Sensitivity:15.71pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:CASP6

Note: This product is for in vitro research use only and is not intended for use in humans or animals.