

Mouse Casp6/Caspase-6 ELISA Kit

Catalog No: #EK20086

Package Size: #EK20086 1x96T



Description

Product Name	Mouse Casp6/Caspase-6 ELISA Kit
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Caspase-6
Target Name	Caspase-6
Other Names	CASP-6,Apoptotic protease Mch-2
Accession No.	O08738

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Casp6

Note: This product is for in vitro research use only and is not intended for use in humans or animals.