

Bovine CRIP1/Cysteine-rich protein 1 ELISA Kit

Catalog No: #EK20096

Package Size: #EK20096 1x96T



Description

Product Name	Bovine CRIP1/Cysteine-rich protein 1 ELISA Kit
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Cysteine-rich protein 1
Target Name	Cysteine-rich protein 1
Other Names	CRP-1
Accession No.	Q56K04

Application Details

Assay Type: SandwichDetect Range:0.625-40ng/mL Sensitivity:0.308ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:CRIP1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.