

Human CRIP1/Cysteine-rich protein 1 ELISA Kit

Catalog No: #EK20097

Package Size: #EK20097 1x96T



Description

Product Name	Human CRIP1/Cysteine-rich protein 1 ELISA Kit
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cysteine-rich protein 1
Target Name	Cysteine-rich protein 1
Other Names	CRP-1,Cysteine-rich heart protein,CRHP,Cysteine-rich intestinal protein,CRIP,CRIP,CRP1
Accession No.	P50238

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.082ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:CRIP1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.