

Mouse Cubn/Cubilin ELISA Kit

Catalog No: #EK20100

Package Size: #EK20100 1x96T



Description

Product Name	Mouse Cubn/Cubilin ELISA Kit
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cubilin
Target Name	Cubilin
Other Names	Intrinsic factor-cobalamin receptor,Ifcr
Accession No.	Q9JLB4

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.11ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Cubn

Note: This product is for in vitro research use only and is not intended for use in humans or animals.