

Human CRAT/Carnitine O-acetyltransferase ELISA Kit

Catalog No: #EK20101

Package Size: #EK20101 1x96T

Description

Product Name	Human CRAT/Carnitine O-acetyltransferase ELISA Kit
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Carnitine O-acetyltransferase
Target Name	Carnitine O-acetyltransferase
Other Names	Carnitine acetylase, Carnitine acetyltransferase, CAT, CAT1
Accession No.	P43155

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.156ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CRAT

Note: This product is for in vitro research use only and is not intended for use in humans or animals.