

Rat Aldh1a7/Aldehyde dehydrogenase, cytosolic 1 ELISA Kit

Catalog No: #EK20102

Package Size: #EK20102 1x96T

Description

Product Name	Rat Aldh1a7/Aldehyde dehydrogenase, cytosolic 1 ELISA Kit
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Aldehyde dehydrogenase, cytosolic 1
Target Name	Aldehyde dehydrogenase, cytosolic 1
Other Names	ALDH class 1,ALDH-E1,ALHDII,Aldehyde dehydrogenase family 1 member A7,Aldehyde dehydrogenase phenobarbital-inducible,Aldh-pb,Aldh1,Aldh1a4
Accession No.	P13601

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.085ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Aldh1a7

Note: This product is for in vitro research use only and is not intended for use in humans or animals.