

Rat Chid1/Chitinase domain-containing protein 1 ELISA Kit

Catalog No: #EK20103

Package Size: #EK20103 1x96T

Description

Product Name	Rat Chid1/Chitinase domain-containing protein 1 ELISA Kit
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Chitinase domain-containing protein 1
Target Name	Chitinase domain-containing protein 1
Accession No.	A0JPQ9

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.082ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Chid1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.