

Human ADRB2/Beta-2 adrenergic receptor ELISA Kit

Catalog No: #EK20105

Package Size: #EK20105 1x96T



Description

Product Name	Human ADRB2/Beta-2 adrenergic receptor ELISA Kit
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Beta-2 adrenergic receptor
Target Name	Beta-2 adrenergic receptor
Other Names	Beta-2 adrenoreceptor,Beta-2 adrenoceptor,ADRB2R,B2AR
Accession No.	P07550

Application Details

Assay Type: SandwichDetect Range:15.6-1000pg/ml Sensitivity:3.9pg/ml Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:ADRB2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.