

Rat HB-beta1 ELISA Kit

Catalog No: #EK20108

Package Size: #EK20108 96T

Description

Product Name	Rat HB-beta1 ELISA Kit
Applications	ELISA
Species Reactivity	Rt

Application Details

Detect Range:0.1 - 16 ng/mL Sample Type: Serum, plasma, Urine,tissue homogenates and other biological fluids.

Sample Volume:50ul

Assay Time: 3-5h

Detection wavelength:450nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.