

Human Melanocortin receptor 4 (MC4R)ELISA Kit

Catalog No: #EK20115

Package Size: #EK20115 96T

Description

Product Name	Human Melanocortin receptor 4 (MC4R)ELISA Kit
Applications	ELISA
Specificity	Homo sapiens human
Other Names	Melanocortin receptor 4,MC4-R,MC4R
Accession No.	UniProtKB:P32245 GenelDoO 4160
Storage	12 months

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.060ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 100 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.