

Chicken D-lactic acid ELISA Kit

Catalog No: #EK20123

Package Size: #EK20123-1 48T #EK20123-2 96T

Description

Product Name	Chicken D-lactic acid ELISA Kit
Applications	ELISA
Species Reactivity	Chicken (Gallus)
Storage	12 months

Application Details

Detect Range:31.25-2000ng/mL

Sensitivity:12.6ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.