

Rat Thrombospondin 1 (THBS1) ELISA Kit

Catalog No: #EK20137

Package Size: #EK20137-1 48T #EK20137-2 96T

Description

Product Name	Rat Thrombospondin 1 (THBS1) ELISA Kit
Applications	ELISA
Species Reactivity	Rt

Application Details

Detect Range:15.6-1000pg/ml

Sensitivity:5.3pg/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.