

Rat Prokineticin 2 (PK2) ELISA Kit

Catalog No: #EK20146

Package Size: #EK20146-1 48T #EK20146-2 96T

Description

Product Name	Rat Prokineticin 2 (PK2) ELISA Kit
Applications	ELISA
Species Reactivity	Rt

Application Details

Detect Range:0.156-10ng/ml

Sensitivity:0.05ng/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.