

# Mouse Sclerostin (SOST) ELISA Kit

Catalog No: #EK20161

Package Size: #EK20161-1 48T #EK20161-2 96T

## Description

|                    |                                   |
|--------------------|-----------------------------------|
| Product Name       | Mouse Sclerostin (SOST) ELISA Kit |
| Applications       | ELISA                             |
| Species Reactivity | Mouse                             |

## Application Details

Detect Range:15.6-1000pg/ml

Sensitivity:6.5pg/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.