

MouseB Beta-CrosslapsB (bCTx)B ELISAB Kit

Catalog No: #EK20169

Package Size: #EK20169-1 48T #EK20169-2 96T



Description

Product Name	MouseB Beta-CrosslapsB (bCTx)B ELISAB Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:78.1-5000pg/ml
Sensitivity:33pg/ml
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.