

General Thyroxine (T4) ELISA Kit

Catalog No: #EK20177

Package Size: #EK20177-1 48T #EK20177-2 96T



Description

Product Name	General Thyroxine (T4) ELISA Kit
Applications	ELISA
Species Reactivity	General

Application Details

Detect Range:1.25-80ng/ml

Sensitivity:0.63ng/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.