

# Human 8-iso prostaglandin (8-iso-PG) ELISA Kit

Catalog No: #EK20178

Package Size: #EK20178-1 48T #EK20178-2 96T



## Description

|                    |                                                |
|--------------------|------------------------------------------------|
| Product Name       | Human 8-iso prostaglandin (8-iso-PG) ELISA Kit |
| Applications       | ELISA                                          |
| Species Reactivity | Hu                                             |

## Application Details

Detect Range:78.1-5000pg/ml

Sensitivity:25pg/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.