

Rat Galactosidase Beta (GLb) ELISA Kit

Catalog No: #EK20190

Package Size: #EK20190-1 48T #EK20190-2 96T

Description

Product Name	Rat Galactosidase Beta (GLb) ELISA Kit
Applications	ELISA
Species Reactivity	Rt

Application Details

Detect Range:31.25-2000pg/mL

Sensitivity:15.2pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.