

# Mouse Aspartate Aminotransferase (AST) ELISA Kit

Catalog No: #EK20200

Package Size: #EK20200-1 48T #EK20200-2 96T



## Description

|                    |                                                  |
|--------------------|--------------------------------------------------|
| Product Name       | Mouse Aspartate Aminotransferase (AST) ELISA Kit |
| Applications       | ELISA                                            |
| Species Reactivity | Mouse                                            |

## Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.118ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.