

Rat Thyroid Stimulating Hormone Receptor Antidoby (TRAb) ELISA Kit

Catalog No: #EK20213

Package Size: #EK20213-1 48T #EK20213-2 96T

Description

Product Name	Rat Thyroid Stimulating Hormone Receptor Antidoby (TRAb) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:125-4000pg/ml

Sensitivity:43pg/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.