

Mouse Lactate Dehydrogenase A (LDHA) ELISA Kit

Catalog No: #EK20234

Package Size: #EK20234-1 48T #EK20234-2 96T



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Mouse Lactate Dehydrogenase A (LDHA) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse
Storage	12 months

Application Details

Detect Range:0.156-10ng/mL
Sensitivity:0.057ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.