

Mouse Defensin-6(DEFA6) ELISA kit

Catalog No: #EK20241

Package Size: #EK20241-1 48T #EK20241-2 96T



Description

Product Name	Mouse Defensin-6(DEFA6) ELISA kit
Applications	ELISA
Species Reactivity	Ms

Application Details

Detect Range:0.312-20ng/mL
Sensitivity:0.126ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.