

Mouse Purinergic Receptor P2Y, G Protein Coupled 12 (P2RY12) ELISA Kit

Catalog No: #EK20254

Package Size: #EK20254-1 48T #EK20254-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Purinergic Receptor P2Y, G Protein Coupled 12 (P2RY12) ELISA Kit
Applications	ELISA
Species Reactivity	Ms
Accession No.	Swiss Prot:Q9CPV9

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.061ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.