

Mouse Reactive OXyen Species (ROS) ELISA Kit

Catalog No: #EK20255

Package Size: #EK20255-1 48T #EK20255-2 96T

Description

Product Name	Mouse Reactive OXyen Species (ROS) ELISA Kit
Applications	ELISA
Species Reactivity	Ms

Application Details

Detect Range:0.156-10U/mL

Sensitivity:0.064U/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.