

Rat Fatty acyl-CoA synthetase (ACS) ELISA Kit

Catalog No: #EK20257

Package Size: #EK20257-1 48T #EK20257-2 96T

Description

Product Name	Rat Fatty acyl-CoA synthetase (ACS) ELISA Kit
Applications	ELISA
Species Reactivity	Rt

Application Details

Detect Range:62.5-4000pg/mL

Sensitivity:25.7pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.