

Mouse Immunoglobulin G2a (IgG2a) ELISA Kit

Catalog No: #EK20258

Package Size: #EK20258-1 48T #EK20258-2 96T



Description

Product Name	Mouse Immunoglobulin G2a (IgG2a) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:0.312-20µg/mL
Sensitivity:0.117µg/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 µL
Assay Time:1-4.5h
Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.