

Rat Glycerol-3-phosphate dehydrogenase, mitochondrial (GPD2) ELISA Kit

Catalog No: #EK20309

Package Size: #EK20309-1 48T #EK20309-2 96T

Description

Product Name	Rat Glycerol-3-phosphate dehydrogenase, mitochondrial (GPD2) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:31.2-2000pg/mL

Sensitivity:12.4pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.