

Rat Interferon Regulatory Factor 7 (IRF7) ELISA Kit

Catalog No: #EK20318

Package Size: #EK20318-1 48T #EK20318-2 96T



Description

Product Name	Rat Interferon Regulatory Factor 7 (IRF7) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:62.5-4000pg/mL

Sensitivity:22.3pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.