

# Sheep Interleukin 8 (IL8) ELISA Kit

Catalog No: #EK20322

Package Size: #EK20322-1 48T #EK20322-2 96T



## Description

|              |                                     |
|--------------|-------------------------------------|
| Product Name | Sheep Interleukin 8 (IL8) ELISA Kit |
| Applications | ELISA Kit                           |
| Storage      | 12 months                           |

## Application Details

|                                                                                                                  |
|------------------------------------------------------------------------------------------------------------------|
| Detect Range:15.6-1000pg/mL                                                                                      |
| Sensitivity:4.21pg/mL                                                                                            |
| Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Sample Volume: 50 uL                                                                                             |
| Assay Time:1-4.5h                                                                                                |
| Detection wavelength:450 nm                                                                                      |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.