

Mouse Hexokinase 2 (HK2) ELISA Kit

Catalog No: #EK20325

Package Size: #EK20325-1 48T #EK20325-2 96T

Description

Product Name	Mouse Hexokinase 2 (HK2) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse
Storage	12 months

Application Details

Detect Range:0.625-40ng/mL

Sensitivity:0.30ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.