

# Mouse Immunoglobulin G (IgG) ELISA Kit

Catalog No: #EK20330

Package Size: #EK20330-1 48T #EK20330-2 96T

## Description

|                    |                                        |
|--------------------|----------------------------------------|
| Product Name       | Mouse Immunoglobulin G (IgG) ELISA Kit |
| Applications       | ELISA                                  |
| Species Reactivity | Mouse                                  |

## Application Details

Detect Range: 1.56-100 µg/mL

Sensitivity: 0.44 µg/mL

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 µL

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.