

Guinea pig Vascular endothelial growth factor A (VEGFA) ELISA Kit

Catalog No: #EK20331

Package Size: #EK20331-1 48T #EK20331-2 96T

Description

Product Name	Guinea pig Vascular endothelial growth factor A (VEGFA) ELISA Kit
Applications	ELISA
Species Reactivity	Guinea pig

Application Details

Detect Range:15.6-1000 pg/mL

Sensitivity:6.3pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.