

# Mouse Rho Associated Coiled Coil Containing Protein Kinase 1 (Rock1) ELISA Kit

Catalog No: #EK20333

Package Size: #EK20333-1 48T #EK20333-2 96T

## Description

|                    |                                                                                |
|--------------------|--------------------------------------------------------------------------------|
| Product Name       | Mouse Rho Associated Coiled Coil Containing Protein Kinase 1 (Rock1) ELISA Kit |
| Applications       | ELISA                                                                          |
| Species Reactivity | Mouse                                                                          |

## Application Details

Detect Range:15.6-1000pg/mL

Sensitivity:5.9pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50  $\mu$ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.