

Mouse Carbohydrate Antigen 15-3 (CA15-3) ELISA Kit

Catalog No: #EK20336

Package Size: #EK20336-1 48T #EK20336-2 96T



Description

Product Name	Mouse Carbohydrate Antigen 15-3 (CA15-3) ELISA Kit
Species Reactivity	Ms

Application Details

Detect Range:3.12-200 mIU/mL
Sensitivity:1.42 mIU/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.